

How Does Microchipping Work?



1. Scan the pet to verify whether he or she has a microchip that has been previously implanted.



2. Scan the Datamars blister pack to verify that the microchip is functional, and that the ID code matches the accompanying barcode labels.



3. Implant the Datamars microchip, then rescan the pet to verify that the microchip can be localized.



4. Do not forget to complete the pet's registration in a database like PetLink!



5. If the pet is ever lost, hopefully his or her microchip can be scanned, and the owner's information retrieved – allowing the pet to find its way home!


creating permanent bonds between owners & pets



PETMAXX™